

## Microlife Device Validation Table

Clinical Devices					
Device	AAMI	BHS	ESH	Circumstance	Reference
Microlife BP 3AS1-2			Pass	At rest	1
Microlife WatchBP O3 (3MZ0)			Pass	At rest	2
Microlife WatchBP O3 (3MZ1)			Pass	At rest	3
Microlife WatchBP Office			Pass	At rest	4
Microlife WatchBP Office ABI		A/A	Pass	At rest	5
			Pass	L-XL size cuff	4,6
				Accuracy in assessing ABI	7
				Accuracy in detecting Atrial Fibrillation	8,9
Microlife WatchBP Office AFIB			Pass	WatchBP Office Equivalence	4*
				Accuracy in detecting Atrial Fibrillation	8, 9
Microlife WatchBP Office Target			Pass	BP A100 Equivalence	10*

Self-measurement Devices					
Device	AAMI	BHS	ESH	Circumstance	Reference
Microlife BP 3AC1-1			Pass	At rest	11
	Pass	A/A		Pregnancy	12
	Pass	B/A		Pre-eclampsia	12
Microlife BP 3AC1-1 PC			Pass	BP 3AC1-1 Equivalence	11*
			Pass	ESRD population	13
Microlife BP 3AC1-2			Pass	BP 3AC1-1 Equivalence	11*
Microlife BP 3AG1		A/A		BP 3BT0-A Equivalence	14*
Microlife BP 3BT0-1		A/A		BP 3BT0-A Equivalence	14*
Microlife BP 3BT0-A		A/A		Small recruitment violation	14
	Pass	A/B		Normotensive Pregnancy	15
	Pass	B/B		Non-proteinuric HBP	15
	Pass	A/B		Pre-Eclampsia	15
Microlife BP 3BT0-A(2)		A/A		BP 3BT0-A Equivalence	15*
Microlife BP 3BT0-AP		A/A		BP 3BT0-A Equivalence	14*
Microlife BP A100			Pass	BP A100 Plus Equivalence	10*



Self-measurement Devices					
Device	AAMI	BHS	ESH	Circumstance	Reference
Microlife BP A100 Plus			Pass	With "M-L Size" Cuff	16
			Pass	At rest	10
			Pass	At rest	17
Microlife Exact BP (BP3MD1-3)			Pass	BP A100 Equivalence M-L cuff	16
Microlife RM 100		A/A		BP 3BT0-A Equivalence	14*
Microlife WatchBP Home (BP3MX1-1)			Pass	At rest	18
	Pass	A/A		Pregnancy	19
	Pass	A/A		Pre-eclampsia	19
Microlife WatchBP Home A (BP3MX1-3)			Pass	Watch BP Home Equivalence	18,19*
				Accuracy in detecting Atrial Fibrillation	8,9
Microlife WatchBP Home N (BP3MX1-4)			Pass	Watch BP Home Equivalence	18,19*
Microlife WatchBP Home S (BP3MX1-5)			Pass	Watch BP Home Equivalence	18,19*
				Accuracy in detecting Atrial Fibrillation	8,9
Microlife WatchBP O3 (3MZ0)			Pass	At rest	2
Microlife WatchBP O3 (3MZ1)			Pass	At rest	3

Wrist Devices					
Device	AAMI	BHS	ESH	Circumstance	Reference
Microlife BP W100	Pass	A/A	Pass	Supine position	20
Microlife BP W200-1			Pass	MAM averaging and supine position	21

\* Indicates that the same validation applies to the device because of equivalence



## References

1. de Greeff A, Nathan H, Stafford N, Liu B, Shennan AH. Development of an accurate oscillometric blood pressure device for low resource settings. *Blood Press Monit* 2008;13:342-8.
2. Stergiou GS, Tzamouranis D, Nasothimiou EG, Karpettas N, Protogerou A. Are there really differences between home and daytime ambulatory blood pressure? Comparison using a novel dual-mode ambulatory and home monitor. *J Hum Hypertens* 2010;24:207-12.
3. Ragazzo F, Saladini F, Palatini P. Validation of the Microlife WatchBP O3 device for clinic, home, and ambulatory blood pressure measurement, according to the International Protocol. *Blood Press Monit* 2010;15:59-62.
4. Stergiou GS, Tzamouranis D, Protogerou A, Nasothimiou E, Kapralos C. Validation of the Microlife Watch BP Office professional device for office blood pressure measurement according to the International protocol. *Blood Press Monit* 2008;13:299-303.
5. Saladini F, Benetti E, Masiero S, Palatini P. Accuracy of Microlife WatchBP Office ABI monitor assessed according to the 2002 European Society of Hypertension protocol and the British Hypertension Society protocol. *Blood Press Monit* 2011;16:258-61.
6. Masiero S, Saladini F, Benetti E, Palatini P. Accuracy of the Microlife large-extra large-sized cuff (32-52 cm) coupled to an automatic oscillometric device. *Blood Press Monit* 2011;16:99-102.
7. Kollias A, Xilomenos A, Protogerou A, Dimakakos E, Stergiou GS. Automated determination of the ankle-brachial index using an oscillometric blood pressure monitor: validation vs. Doppler measurement and cardiovascular risk factor profile. *Hypertens Res* 2011;34:825-30.
8. Wiesel J, Fitzig L, Herschman Y, Messineo FC. Detection of atrial fibrillation using a modified microlife blood pressure monitor. *Am J Hypertens* 2009;22:848-52.
9. Stergiou GS, Karpettas N, Protogerou A, Nasothimiou EG, Kyriakidis M. Diagnostic accuracy of a home blood pressure monitor to detect atrial fibrillation. *J Hum Hypertens* 2009;23:654-8.
10. Stergiou GS, Giovas PP, Neofytou MS, Adamopoulos DN. Validation of the Microlife BPA100 Plus device for self-home blood pressure measurement according to the International Protocol. *Blood Press Monit* 2006;11:157-60.
11. Topouchian JA, El Assaad MA, Orobinskaia LV, El Feghali RN, Asmar RG. Validation of two devices for self-measurement of brachial blood pressure according to the International Protocol of the European Society of Hypertension: the SEINEX SE-9400 and the Microlife BP 3AC1-1. *Blood Press Monit* 2005;10:325-31.
12. de Greeff A, Reggiori F, Anthony J, Shennan A. The Microlife 3AC1: An accurate blood pressure measurement device in pregnancy and pre-eclampsia. *J Hypertens* 2006;24.
13. Thompson AM, Eguchi K, Reznik ME, Shah SS, Pickering TG. Validation of an oscillometric home blood pressure monitor in an end-stage renal disease population and the effect of arterial stiffness on its accuracy. *Blood Press Monit* 2007;12:227-32.
14. Cuckson AC, Reinders A, Shabeeh H, Shennan AH. Validation of the Microlife BP 3BTO-A oscillometric blood pressure monitoring device according to a modified British Hypertension Society protocol. *Blood Press Monit* 2002;7:319-24.
15. Reinders A, Cuckson AC, Lee JT, Shennan AH. An accurate automated blood pressure device for use in pregnancy and pre-eclampsia: the Microlife 3BTO-A. *BJOG* 2005;112:915-20.
16. Bonso E, Dorigatti F, Palatini P. Accuracy of the BP A100 blood pressure measuring device coupled with a single cuff with standard-size bladder over a wide range of arm circumferences. *Blood Press Monit* 2009;14:216-9.
17. Belghazi J, El Feghali RN, Moussalem T, Rejdych M, Asmar RG. Validation of four automatic devices for self-measurement of blood pressure according to the International Protocol of the European Society of Hypertension. *Vasc Health Risk Manag* 2007;3:389-400.
18. Stergiou GS, Giovas PP, Gkinos CP, Patouras JD. Validation of the Microlife WatchBP Home device for self home blood pressure measurement according to the International Protocol. *Blood Press Monit* 2007;12:185-8.
19. Chung Y, de Greeff A, Shennan A. Validation and compliance of a home monitoring device in pregnancy: microlife WatchBP home. *Hypertens Pregnancy* 2009;28:348-59.
20. Palatini P, Dorigatti F, Bonso E, Ragazzo F. Validation of Microlife BP W100 wrist device assessed according to the European Society of Hypertension and the British Hypertension Society protocols. *Blood Press Monit* 2009;14:41-4.
21. Palatini P, Dorigatti F, Bonso E, Ragazzo F. Validation of the Microlife BP W200-1 wrist device for blood pressure measurement. *Blood Press Monit* 2008;13:295-8.